

## MILESTONE 1 (STAGE 2) – LIST OF OBJECTIVES, CONSTRAINTS, AND FUNCTIONS

Team ID: Mon-16

1. As a team, create a final a list of objectives, constraints, and functions in the table below.

- Use your individual *Pre-Project Assignment* to build your team's final list
- The exact number you should have depends on what information you have gathered from the Project Module.

Objectives	Constraints	Functions
Container should be rigid and durable	<i>Container must fit within the assigned autoclave footprint</i>	Securely holds surgical tools
Resistant to high temperature	<i>Base must lay parallel to the assigned autoclave floor</i>	Facilitates sterilization
Unreactive with chemicals used to sterilize	<i>Must be able to be picked up by the robot arm's end effector – max width of 80-150 mm</i>	Secures tools during travel
Low thermal expansion	<i>All features must exceed 2 mm</i>	Be able to be picked up by robot arm
Identifiable by the robot arm	<i>Mass does not exceed 350g</i>	Store liquid
Impervious to liquids		

2. What is the primary function of the entire system?

Allows the sterilization of surgical tools.
---

3. What are the secondary functions?

Securely holds surgical tools
Able to be picked up by robot arm
Store liquid

## MILESTONE 1 (STAGE 3) – MORPHOLOGICAL ANALYSIS

Team ID: Mon-16

1. Identify multiple means to perform the secondary functions that your team came up with during Stage 1 of this milestone. One sub-function (pick up) is already listed for you. The other two sub-functions are for your team to choose.

→ Make sure that every mean for the “pick up” sub-function assumes that the end effector of the robot arm is a gripper. The means for your other sub-functions do not need to follow this assumption.

Function	Means					
Pick up	Handle	Holes for gripper	Hook	Lip on the edge		
Tool Storage	Box	Basket	Tool cutout	Adhesive strips	Magnet	Pouch
Store liquid	Bowl	Container	Cup	Box		

## MILESTONE TWO (TEAM): SUBTEAMS, SKETCHES, & WORKFLOW

### PROJECT TWO: MILESTONE 2 – COVER PAGE

Team ID: Mon-16

Please list full names and MacID's of all *present* Team Members

Full Name:	MacID:
Mohammad Bilal	bilalm14
Owen Brazel	brazelo
ZiDi Yao	yaoz25
Ibrahim Arif Qadeer	qadeeri
Muhammad Saad Siddiqi (Online)	siddim98

## MILESTONE 2 (STAGE 3A) – LOW-FIDELITY PROTOTYPE (MODELLING SUB-TEAM)

Team ID: 

Mon-16
--------

Complete this worksheet during design studio 8 after creating the low-fidelity prototypes.

1. Take multiple photos of your low-fidelity prototypes  
→ Include an index card (or similar) next to the prototype, clearly indicating your Team Number, Name and MacID on each sketch
2. Insert your photo(s) as a Picture (Insert > Picture > This Device)
3. **Do not include more than two prototype photo's per page**

Make sure to include photos of <u>each</u> team member's prototype
--

Team ID: Mon-16

Name: Muhammad Saad Siddiqi

MacID: siddim98

*Insert screenshot(s) of your low-fidelity prototype below*



*Insert screenshot(s) of your low-fidelity prototype below*





Team ID: 

Mon-16

Name: Owen Brazel


MacID: brazelo

*Insert screenshot(s) of your low-fidelity prototype below*





## ENGINEER 1P13 – Project Two: *Get a Grip*



\*If you are in a sub-team of 3, please copy and paste the above on a new page



## MILESTONE 2 (STAGE 3B) – LOW-FIDELITY PROTOTYPE OBSERVATIONS (MODELLING SUB-TEAM)

Team ID: Mon-16

As a sub-team, document your observations for each low-fidelity prototype. Make sure to label your observations to indicate which prototype it belongs to. As a starting, consider the following: (note, this does not fully encompass all discussion points)

- Advantages and disadvantages of each prototype
- Extent to which each concept aligns (or does not align) with the List of Objectives, Constraints, and Functions you came up with for Milestone 1
- Reliability of the design in picking up the surgical tool
- Reliability of the design in securing the surgical tool
- Extent to which it allows for tool sterilization

*Document your observations for each prototype in the space below. It is recommended you document observations in a **table** or in bullet form (it should be clear which prototype you are referring to for each observation).*

	Prototype 1 – Designed by Owen Brazel (Low-fidelity prototype and observations made by Saad Siddiqi)	Prototype 2 – Designed by Saad Siddiqi (Low-fidelity prototype and observations made by Owen Brazel)
Advantages	<ul style="list-style-type: none"> <li>- Holes around the container allow for sterilization to occur smoothly</li> <li>- Cut out at the bottom securely holds the surgical tool</li> <li>- Lip on outer edge aids the robot arm to transfer the container easily</li> </ul>	<ul style="list-style-type: none"> <li>- Allows for easy sterilization</li> <li>- Securely holds the tool</li> <li>- Can hold tool easily</li> <li>- Has a lip on the outer edge of the box to allow the robot arm to pick it up easier</li> </ul>
Disadvantages	<ul style="list-style-type: none"> <li>- No upper lid or strap to stop tool from falling out of the cut out if the container is flipped over</li> </ul>	<ul style="list-style-type: none"> <li>- No handle</li> <li>- No lid to stop the tool from falling if the tool falls from the holder</li> </ul>
Objectives (Not Aligned)	<ul style="list-style-type: none"> <li>- Impervious to liquids</li> </ul>	<ul style="list-style-type: none"> <li>- Impervious to liquids</li> </ul>
Constraints (Aligned)	<ul style="list-style-type: none"> <li>- Max width of 80-150 mm</li> <li>- All features exceed 2 mm</li> <li>- Mass does not exceed 350 g</li> </ul>	<ul style="list-style-type: none"> <li>- Edges are more than 2mm</li> <li>- Design allows for transportation</li> <li>- Box is not more than 350g</li> <li>- Max width of 80-150mm</li> </ul>
Functions (Not Aligned)	<ul style="list-style-type: none"> <li>- Store liquid</li> </ul>	<ul style="list-style-type: none"> <li>- Store liquid</li> </ul>

## ENGINEER 1P13 – Project Two: *Get a Grip*

Reliability	<ul style="list-style-type: none"><li>- Design is reliable in picking up surgical tool</li><li>- Design is reliable in securing surgical tool</li><li>- Design allows for tool sterilization</li></ul>	<ul style="list-style-type: none"><li>- Design is reliable in picking up surgical tool</li><li>- Design is reliable in securing surgical tool</li><li>- Design allows for tool sterilization</li></ul>
-------------	--	--

## MILESTONE 2 (STAGE 4A) – WORKFLOW PEER-REVIEW (COMPUTATION SUB-TEAM)

Team ID: Mon-16

As a sub-team, document your observations, specifically any similarities and differences between each team member's visual storyboard or flowchart in the table below.

*Document your observations for each visual storyboard / flowchart in the space below.*

Similarities:

- Both flowcharts start off at home position and move to pick up location
- Both flowcharts involve a decision to be made to determine the size of the container
- Both flowcharts perform the same function based on the size of container (open and close drawer if the container is large)
- Both Flowcharts involve looping if all containers have not been moved

Differences:

- One flowchart uses potentiometers to adjust the grippers
- One flowchart starts off with a decision whether there are any containers left
- One flowchart does not check at the end whether there are any other containers left to move

## MILESTONE 2 (STAGE 4B) – PROGRAM PSEUDOCODE (COMPUTATION SUB-TEAM)

Team ID: Mon-16

As a sub-team, write out a pseudocode outlining the high-level workflow of your computer program in the space below.

Program starts

While all containers have not been moved:

    Q-arm moves to home position

    Move Q-arm to pickup location using predetermined xyz coordinates

    Using the right potentiometer to control the angle of rotation.

    Using the left-potentiometer to control which position to place in the autoclave

    If container size matches the position of the autoclave(determine through container id):

        Then move q-arm to drop off location using xyz coordinates

        Open autoclave drawer

        Move q-arm back to container location

        Pick up container and move to drop off location

        Close autoclave drawer

    Elif container is small (determine through container id):

        Then pickup container and move to drop off location

    If all containers have been moved:

        Break out of loop

Program ends